

Imidazole Buffered Saline

CE

Intended Use

Imidazole Buffered Saline is a buffer that is intended to be used to dilute plasma samples for fibrinogen determination or other testing.

REAGENT

FOR IN-VITRO DIAGNOSTIC USE ONLY

Ingredients: Each vial contains 135 mL of imidazole Buffer. Imidazole Buffer contains 15mM Imidazole, 0.125 M Sodium Chloride and 0.02% Sodium Azide as a preservative.

WARNING: Sodium Azide: This buffer is preserved with sodium azide which can form highly explosive metal azides if exposed to lead or copper in plumbing. This buffer should be discarded into a sink only with large volumes of water to minimize such risk.

Preparation for Use: The buffer is packaged ready for use.

Storage and Stability: The buffer is stable until the expiration date printed on the label when stored at 2-8 ° C.

LABORATORY TESTING

If the Imidazole Buffer is to be used as a diluent for a fibrinogen assay, it should be used in accordance with the laboratory procedures for fibrinogen determination that are in place for each individual laboratory.

